

# Sample Collection Kit Checklist

<b>Urine Specimen Collection Kit</b>	<b>QTY</b>	
Urine cup	1	_____
Biohazard ziploc bag	1	_____
Metals-free tube w/ sulfamic acid, 15mL ( <i>Hg aliquoting</i> )	1	_____
Syringe, 10 mL slip tip ( <i>Hg aliquoting</i> )	1	_____
Nitrile powder-free gloves	2	_____
Divider containers	1	_____
Adsorbent pad	1	_____
Rack for falcon tubes	1	_____
Cooler	1	_____
Sealing tape		_____
Dry ice		_____
Labels		_____
Towel		_____

<b>Water Sample Collection Kit</b>	<b>QTY</b>	
250mL HDPP container	1	_____
As speciation cartridge	1	_____
Syringe, 30 mL slip tip ( <i>As speciation</i> )	1	_____
250mL HDPP container w/ Aluminum foil ( <i>As speciation</i> )	1	_____
Nitrile powder-free gloves	2	_____
Cooler	1	_____
Dry ice/Ice pack		_____
Sealing tape		_____
Labels		_____

<b>Administrative</b>	<b>QTY</b>	
Patient Consent Form	1	_____
Participant Collection Checklist	1	_____
Participant Collection Instructions	1	_____
Urine Sample Information Sheet	1	_____
Sample Collection Process Checklist	1	_____